SHAPE YOUR PERFECT BODY




LET THE COLD TRANSFORM YOUR LIFE





¿HOW DOES cooltech WORK?






TWO SIMULTANEOUS APPLICATORS

The exclusive design of **cool**tech means that you can work at a temperature of **with two applicators simultaneously**, thereby reducing treatment time and favouring the patient's comfort, since in just one session you will be able to reduce two areas, obtaining a more visible effect.

The **multiple body areas** that can be treated with **cool**tech and its special applicators is another one of the enormous advantages of this equipment, since full body contouring can be done in just a few sessions.





Cool Gel Pad

PATENTED MEMBRANE

cocoon medical has patented the cool gel pad membrane, which guarantees maximum safety for its patients

Its mandatory use allows carrying out treatments with the best results and the greatest skin protection.

#safetyandresults







cooltech APPLICATORS AND COUPLERS

To improve the access of **cool**tech applicators to reduced fat folds or hard fat, we have developed our own couplers that facilitate treatment as well as allowing access to more areas and for more clients.





Treatments indicated for Straight HP

It is the head which has greatest capacity. Indicated for treating large areas on flat surfaces. (Folds thicker than 35 mm and lengths up to 15 cm).





Treatments indicated for Curved HP

Specially designed for working on curved areas without losing suction capacity. (Folds thinner than 35 mm and lengths up to 15 cm).





Treatments indicated for Double HP

Designed to be able to treat the entire length of the abdominal area in a single session and in a homogeneous way. (Folds thicker than 35 mm and lengths up to 21 cm).





Treatments indicated for Tight HP

Designed for treating narrow fat folds. This is the ideal complement for the Straight HP in final sessions, being able to completely remove more localised fat. (Folds thinner than 35 mm and lengths up to 15 cm).





Treatments indicated for Tiny HP

Its small size is indicated for areas of highly reduced and difficult to access fat, such as axillary fold, double chin, upper fold of the knees, etc. (Folds thinner than 20 mm and lengths up to 7 cm).





Treatments indicated for Oval HP

Its wide head opening allows treating large areas and hardened fat in which suction is difficult. Ideal for legs and flanks. (Areas of $10 \times 15 \text{ cm}$).



DOCTORS' TESTIMONIES

cooltech is by far my favourite method for reducing fat, because its techniques are highly customisable and it has different applicators for each area, that way being able to obtain truly miraculous results. My clients and I are truly delighted with this piece of equipment."

Dra. Lohia (North American dermatologist specialized in treating celebrities)

cooltech is the perfect method for reducing localised and general fat in patients who do not wish to undergo surgical procedures. Its variety of applicators mean that all kinds of patients can be treated, always with incredible results." Dr. Serena (Cosmetic doctor with more than 20 years' experience in the sector)

SATISFIED PATIENTS AND 95% OF OBJECTIVE RESULTS

PATIENTS' TESTIMONIES

I am very happy, a month after the treatment I was able to buckle my trousers without my abdominal fat hanging over."

J. Ford. 56. Abdominal zone

I am satisfied with the treatment. After two sessions, I can still feel how it has reduced the fat I had on my abdomen. I like the results and I think it is recommendable."

C. Matías. 49. Abdominal zone



Before

DATE: 05/07/16 CONTOUR: 90 cm



After: 2 sesiones

DATE: 15/12/16 CONTOUR: 82 cm



Before DATE: 28/09/16 CONTOUR: 103 cm



After: 1 sesión **DATE:** 18/11/16



CONTOUR: 96 cm




Before

DATE: 05/09/16 CONTOUR: Izquierda: 66 cm Derecha: 66 cm



After: 1 sesión

DATE: 10/11/16 CONTOUR: Izquierda: 61 cm Derecha: 62 cm





Before DATE: 04/07/16 CONTOUR: 99.5 cm



After: 2 sesiones **DATE:** 07/11/16 CONTOUR: 93 cm





cooltech

by cocoon medical

cooltech is a new concept in cryoadipolysis treatment. Its various applicators cool the area thereby achieving extraordinary results from the very first session.



Double **HP**

Oval HP

www.cocoonmedical.com











